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about us:

The Polish Geological Institute is a research institute having the status of a National Research Institute. It was established in May 1919 and is supervised by the Minister of the Environment. Under the Act on Geological and Mining Law, and the Act on Water Law, it performs the tasks of the Polish Geological Survey and the Polish Hydrogeological Survey.

The Polish Geological Institute – NRI performs its mission through intense activities in all fields of earth sciences all over Poland. It is the main depositary and a source of knowledge, information and geological, hydrogeological and geoenvironmental data in Poland. The Association of the Geological Surveys of Europe, EuroGeoSurveys, will remain the basic platform of international cooperation for the Institute. All samples prior to performing measurements on SHRIMP microprobe shall be fully documented photographically using four types of images:


- reflecting light optical imaging () 1)
- transmitted light optical imaging () 2)
- SEM –CL monochromatic cathodoluminescence imaging or BSE imaging () 3)
- low magnification -SEM imaging to assembly a navigation map (4)

### SAMPLE PREPARATION

These images allow the complexity of the internal structure of the sample to be revealed, a type of material, the quality of its surface, the existence of the inclusions and consequently precisely select the area for SRIMP measurements.

**~**α

### SHRIMP Ile/MC | SAMPLE PREPARATION


Isotopic studies are carried out on samples such as minerals, fossils and other items of very small size. They are embedded in an epoxy resin together with the proper standards, and then are cut to expose the surface for analysis and polished precisely. This creates a sample mount.

Monomineral samples for ion microprobe SHRIMP IIe/MC measurements can be prepared in the sample separation facility in PGI-NRI. It is provided with the necessary equipment to complete crushing and mineral separation, which allows the proper sample preparation. The SHRIMP accepts two styles of epoxy mounts: A standard epoxy resin of 25.4 mm () a) diameter and a "mega-mount" that is slightly larger at about 35 mm () b). This stage of sample preparation should be done at PGI. There is also the possibility of working on a standard thin sections () c).

#### **APPLICATION**

#### TRACE ELEMENTS

SHRIMP IIe/MC ion microprobe -the ion mode and elements detection

Oxygen ( O- 2) primary source: Positive secondary ion mode Caesium (Cs+) primary source: Negative secondary ion mode Quantifiable in either positive and negative ion mode

shrimp@pgi.gov.pl

TRACE ELEMENTS

SHRIMP analysis has been applied to numerous problems in trace element geochemistry, such as:

- measurements of the concentration and distribution of trace elements within individual grains, and mineral solid inclusions
- substitution mechanisms within individual grains
- rare elements in geological and other samples
- depth profiling of surface films and diffusion profiles and isotopic diffusion rates
- inclusions wthin steel
- partitioning of elements between phases, and dissolution rates of minerals
- isotope imaging by rastered sample

### APPLICATION MULTI -ISOTOPE ANOMALIES

**Baszkówka meteorite - L5** ordinary chondrite, weighing of 15.7 kg, and dropped at 25.VIII.1994 in Poland, 52 ° 02'N, 20 ° 56'E. Chondrites were formed by the accretion of particles of dust and grit present in the primitive Solar System which gave rise to asteroids.

**Baszkówka meteorite** is oriented meteorite, with perfectly preserved distinctive shape of a flattened cone with radiating grooves. The entire surface is covered by a black melted crust.

MULTI -ISOTOPE ANOMALIES

The SHRIMP IIe/MC can be used to determine isotopic ratios of elements on analytical spots of similar size.

Although the SHRIMP lle can achieve sub-permil level precision using peak-switching analyses under ideal conditions, use of the multicollector is advised for maximizing performance.

- solar wind isotope analysis of lunar materials
- examining stellar nucleosynthesis
- investigating Ti isotopic ratios in meteorites
- determining Pb isotopic composition of lunar samples
- anomalies in refractory meteoritic minerals
- the analysis of radionuclide material (U and Pu isotope ratios) for nuclear forensics
- anomalies in experimental samples

#### The SHRIMP ion microprobe is used for stable isotope determination in biogenic and inorganic minerals for:

- conodont thermometry on biogenic apatite as a record of the sea water paleotemperature
- tracing climate change "subtle paleoclimatology" on geological timescales on CO<sub>3</sub> bearing fossils
- the impact of changes in sea water paleotemperature on biodiversity
- the exploration of isotopic heterogeneities of mammalian tooth enamel for investigating the diets of fossil taxa

#### **APPLICATION**

STABLE ISOTOPES

- analysis of the stable isotope ratios of teeth and of bone of ancient people as an indicators of palaeoclimate, palaeodiet, and palaeoecology
- seasonal oxygen isotope cyclicity and annual time markers of osteogenesis from polished bone cross-sections
- sulphur isotopes genetic study on micro scale in the sulphide minerals that form massive sulphide deposits
- carbon isotopic composition
- lithium and boron isotopic changes on clay minerals related to hydrocarbon migration


## SHRIMP IIe/MC APPLICATION

#### STABLE ISOTOPES

The new advances of SHRIMP ion microprobe have been in the field of light stable isotope ratios, especially oxygen, sulphur and carbon. Stable isotope geochemistry examines the change of the isotopic composition of an elements (H, Li, B, C, N, O, S) produced by chemical or physical processes. The stable isotope compositions of biogenic materials record a combination of the biological processes and environmental parameters. Oxygen isotopic compositions of individual spots can be measured in less than 5 minutes with internal precisions of about 0.1 per mille. Traditionally, oxygen isotopic compositions are measured following chemical processing of large samples weighing several milligrams. The SHRIMP IIe/MC ion microprobe allows exploration of isotopic heterogeneities on a very small scale, for example the composition of individual zoning in microfossils.

Moreover SHRIMP measurement of S isotopes has been used to understand mineral growth mechanisms, follow changes in fluid compositions and to constrain the conditions under which host rocks and ore deposits form.

conodonts

### APPLICATION

ADVANCED GEOCHRONOLOGY

1 mm

The zircon grains monomineral concentrate – a basic material for U-Pb SHRIMP isotopic dating SHRIMP IIe/MC APPLICATION

#### ADVANCED GEOCHRONOLOGY

Through U-Pb dating of zircons and others accessory U-Pb-Th bearing minerals such as zircon, titanite, perovskite, allanite, rutile, baddeleyite, monazite, and xenotime and opal where U-Th disequilibrium geochronology has also been performed.

SHRIMP analyses are advantageous in that the fine spatial resolution allows targeting of subgrain domains in mineral grains with complex growth histories, avoiding U or Pb-rich inclusions or metamict domains. When temporal precision higher than 1% is required, the SHRIMP IIe/MC can be used to screen zircon grains for lead loss, inherited cores, or other problems that may inhibit or complicate TIMS analysis.

The small analytical volume (~1000 cubic microns) allows minimal sample consumption prior to subsequent analytical techniques. Therefore SHRIMP IIe is a appropriate tool for:

- accurate determination of the ages of igneous rocks
- accurate determination of the ages of mineral deposits
- magma provenance from inherited zircon cores and age of metamorphism from zircon overgrowths
- unravelling the history of complex metamorphic terranes,
- tracing crustal growth and recycling through geologic time
- sediment provenance and correlation using detrital zircon cores
- deposition ages from diagenetic overgrowths
- tectonic framework and timing of basin formation for petroleum exploration
- calibrating the Palaeozoic time-scale
- dating of the Earth's oldest crust
- examining the oldest zircons in the solar system

## SHRIMP Ile/MC APPLICATION


GEOLOGY, CHEMISTRY, COSMOCHEMISTRY, (PALAEO-) BIOLOGY, NUCLEAR ENERGY, PHYSICS, MATERIALS SCIENCES,

SHRIMP ion microprobe originally was designed for geological research purposes. It was mainly focused on geochronology. The first instrument was constructed by the Research School of Earth Sciences (RSES) at the Australian National University (ANU) in Canberra and then SHRIMP was developed and manufactured by Australian Scientific Instruments (ASI).

The application of SHRIMP IIe/MC, however, is not limited to the precise U-Pb isotopic geochronology. A relatively new application of SHRIMP IIe is measurement of the stable isotopes of the light elements including: oxygen (O), sulfur (S) and carbon (C). Without chemical dissolution and destruction, in-situ, microanalysis of nano amount of sample for diagnosis (only 2 to 5 nanograms) the SHRIMP IIe/MC enables detection of the isotopic variations within different specimens.

The SHRIMP ion microprobe is a versatile and most efficient microanalytical tool for a wide range of geochemical and isotopic investigations. SHRIMP can be used to solve crucial problems for geological, environmental, or technological and experimental studies.

For many years SHRIMP has provided the world's geological surveys and the mining industry with valuable isotopic data to reduce exploration risk. There are two main types of analysis which dominate current SHRIMP application, which are: U-Pb-Th geochronology and stable isotopes ratios of the light elements.

The quality of the instrument can be assessed by its ability to detect trace elements present in the target at low concentrations (sensitivity) and its ability to distinguish between ions of very similar mass (mass resolution).

The high mass resolution of SHRIMP IIe is achieved by the use of a double-focusing mass spectrometer (energy and mass refocusing) with a very large turning radius of magnet and electrostatic analyzer (magnet radius 1000 mm, electrostatic analyzer radius 1272 mm).


Therefore SHRIMP has a beam line over 7 m long and weighs more than 13 tons.

### TECHNICAL PARAMETERS

The high sensitivity at high mass resolution:  $5400 \text{ M}/\Delta M$  at 80% transmission with flat – topped peaks, permits resolution of major molecular interferences during analysis.

SHRIMP can be used for a routine measurement of the isotopic composition of the abundances of most elements in the Periodic Table and can operate in positive or negative ion mode.

SHRIMP is equipped with a Duoplasmatron for positive secondary ion analysis and a Caesium source and charge neutralization system for negative secondary ion analyses.

Caesium Gun, which provides a stable positive Cs+ beam

 View of the opened source chamber of the SHRIMP IIe/MC, with the primary column and Köhler apertures (location showing by hand) and secondary column extending (to the left) towards the electrostatic analyser


#### shrimp@pgi.gov.pl

### SHRIMP Ile/MC TECHNICAL PARAMETERS

SHRIMP makes in situ isotopic and chemical analyses of complex solid materials (without chemical dissolution) by bombarding the sample surface with a high energy primary ion beam in a high vacuum, with a spot diameter of only a few microns.

This desorbs surface species through a physical process called sputtering. The primary ion beam diameter can be set from 30 microns to less than 5 microns. The sputtered fragments are ionised, and then the secondary ions are gathered using electrostatic lenses and transferred to a mass spectrometer in which they are separated according to their relative masses.

The secondary mass analyser comprises an electrostatic sector, a quadrupole correcting lens and electromagnet. This design gives small values for all second-order image aberrations.


Interchangeable sources of the primary ions:

Duoplasmatron, which produce negative ion beam from pure O2


SHRIMP IIe/MC is an extremely sensitive tool for micro-scale chemical and isotopic analysis which allows us to take a look back billions of years in time and respond to questions concerning the geological history of the Earth and early stages formation of the crust, or even confront the history of the solar system.

Ministry of Science and Higher Education Republic of Poland

The **S**ensitive **H**igh **R**esolution Ion **M**icro**P**robe (SHRIMP) is a secondary ion mass spectrometer (SIMS) with high sensitivity at high mass resolution. SHRIMP utilizes a beam of charged particles to probe solids (such as individual mineral grains), for the purpose of elemental and isotopic analysis. SHRIMP ion microprobe as a significant instrumental investment in PGI-NRI is made possible by a grant from the Fund for Polish Science and Technology (661/FNiTP/616/2011) founded by the Polish Ministry of Science and Higher Education.

### Ion Microprobe Facility SHRIMP IIe/MC

The Ion Microprobe Facility is part of the Micro-area Analysis Laboratory, which is dedicated to qualitative and quantitative analysis of the chemical composition on solid material using micro- amount of samples.

The original idea to build that highly specialized facility in the Institute was put forward by Director Prof. Jerzy Nawrocki. Thanks to grant from the Polish Science and Technology Fund of the Minister of Science and Higer Education, it became possible to equip the facility in the latest, upgraded version of SHRIMP IIe/MC sensitive and high-resolution ion microprobe with multicollector.

This instrument has been manufactured by Australian Scientific Instruments (ASI), with cooperation of the Australian National University in Canberra, where the ion microprobe was first designed and developed. The SHRIMP instrument has allowed the discovery the oldest minerals on Earth - zircons from Jack Hills, dated at 4.4 billion years. The PGI ion microprobe is the fourth of the SHRIMP instruments in the world designed to a new enhanced configuration.

Resarche works in the new facility will focus mostly on geochronology and determination of the isotopic ages of rocks and minerals, as crucial for many fields of geology. The technical parameters of this unique and versatile device allow however, for much more. Currently the high sensitivity stable isotope analysis enables the determination of climate, changes during the last stages of the Earth's history. SHRIMP'S perfect technology is still evolving and will certainly soon provide a much surprising data for various fields of modern geology.

#### contact:

#### Micro-area Analysis Laboratory Ion Microprobe Facility

#### Polish Geological Institute National Research Institute

shrimp@pgi.gov.pl

4, Rakowiecka Street 00-975 Warsaw, Poland www**.pgi**.gov.pl

### MICRO-AREA ANALYSIS LABORATORY

### SHRIMP IIe/MC ION MICROPROBE SECTION

Polish Geological Institute National Research Institute Polish Geological Survey Polish Hydrogeological Survey

shrimp@pgi.gov.pl